

C-Tag™ Technology //// Fast and Accurate Results  
Multiplex Structure //// Semi-Quantitative Detection  
High Sensitivity and Specificity /// Early Diagnosis  
Automatic and Easy Analysis with SCM-PLEX™ Viewer  
Multiplex Structure //// Semi-Quantitative Detection  
Automatic and Easy Analysis with SCM-PLEX™ Viewer  
SCM-PLEX™ C-Tag™ Technology / Early Diagnosis  
C-Tag™ Technology //// Fast and Accurate Results  
Multiplex Structure //// Semi-Quantitative Detection  
High Sensitivity and Specificity /// Early Diagnosis



info@selcouthlifesciences.com •  
M: +91 85117 92850



**SELCOUTH**  
Lifesciences

Powered by



SCM Plex™  
**HPV29**  
Detection Kit



Easy and Fast

Easy and Fast



# SCM Plex™ HPV29 Detection

## Simultaneous Genotyping of 29 HPV Genotypes

High-risk HPV	16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 68, 73, 82 (15 types)
Probable high-risk HPV	26, 53, 66, 67, 69 (5 types)
Low-risk HPV	6, 11, 40, 42, 43, 44, 54, 61, 70 (9 types)

- **Single-reaction Real-time PCR Assay**  
By SCM Biotech's proprietary C-Tag™ technology
- **Superior Clinical Performance**  
Clinical Sensitivity: 97.53% - 100%  
Clinical Specificity: 99.81% - 100%  
Reproducibility: 100%
- **Ease of Use**  
Automated analysis with SCM Plex™ Viewer of software  
Whole process monitoring with internal control  
UDG system to prevent carry-over contamination
- **Fast Turnaround Time (TAT)**  
Within 2.5 hours from nucleic acid extraction to result analysis,  
Same-Day Report

